

Appendix: 3 easy steps to convert counts per minute (CPM) to Becquerel per centimeter squared (Bq/cm²)

Step 1: Determine the Radionuclide Class

The CNSC has grouped radionuclide into [3 classes](#). The 3 classes are shown below:

Class A Radionuclide:

Cs-134, Cs-137, Na-22, Po-210, U-235, U-238

Class B Radionuclide: are long-lived or emit beta or gamma radiation

Co-58, Sr- 90, Y-90

Class C Radionuclide: are short-lived and emit beta or gamma radiation

C-14, Ca-45, Cl-36, Cr-51, Fe-55, H3, I-125, Ni-63, P32, P-33, S-35, Tc-99, Tc-99m

Step 2: Determine the regulatory limit

Each of these has limits assigned in terms of unfixed contamination in either a radioactive use or storage area, or in an area where radioactive material is used.

Non - Fixed Contamination	Class A radionuclide	Class B radionuclide	Class C radionuclide
in all areas, rooms or enclosures where unsealed nuclear substances are used, or stored	3 Bq/cm ²	30 Bq/cm ²	300 Bq/cm ²
in all other areas and packaging	0.3 Bq/cm ²	3 Bq/cm ²	30 Bq/cm ²

Step 3: Convert CPM to Bq/cm²

The readings from contamination meters and non-portable instruments are related to regulatory criteria if the efficiency of the instrument for a specific radioisotope is known. Instrument efficiencies for specific radioisotopes can be obtained from the manufacturer or determined using an appropriate standard of known activity.

How is 'Becquerel per centimeter squared' calculated when using a Liquid Scintillation Counter/Gamma Counter (wipe testing)?

$$\text{Result in Bq/cm}^2 = (\text{Reading in cpm} - \text{bkg}) / (E_c \times E_w \times 60 \times A)$$

where bkg = counts per minute of the background,
E_c = counter efficiency (refer to vendor's manual),
E_w = wipe efficiency, assume 10% (0.1), and
A = area wiped in cm².

How is 'Becquerel per centimeter squared' calculated when from CPM when using a survey meter?

$$\text{Result in Bq/cm}^2 = (\text{Reading cpm} - \text{bkg}) / (E_c \times 60 \times A)$$

where bkg = counts per minute of the background,
E_c = survey meter efficiency (refer to vendor's manual), and
A = area wiped in cm² (19.6 cm² for a pancake probe).