




Fixed-Angle Rotor, Aluminum

Major applications: Harvesting bacteria, processing tissue homogenates, subcellular particulates, routine pelleting such as precipitates and phase separations.

	Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
Outer Row	25000	74200	265	24 x 15 mL 18 x 100 mm	360 mL
Inner Row	25000	60200	380	0.75 x 4 in	

For use in Avanti J Series centrifuges
(except Avanti J-E) centrifuges.

No. 363050. JA-25.15 Fixed-Angle Rotor with Biosafety Lid.
No. 363054. JA-25.15 Fixed-Angle Rotor, Single Lock Lid.

Tubes and Bottles

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Cavity Diameter (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed (Outer Row)
Bottles with Screw-On Caps									
Polycarbonate	10.0	355672	25	16 x 80	342327	1	60 200	204	25 000
Open-Top Tube									
Thickwall Polypropylene	4.0	355644	25	13 x 64	303313	1	49 000	120	25 000
	10.0	355640	25	16 x 75	342327	1	56 700	180	25 000
Thickwall Polycarbonate*	4.0	355645	25	13 x 64	303313	1	49 000	120	25 000
	10.0	355630	25	16 x 75	342327	1	56 700	180	25 000
	15.0	342080	100	18 x 98	—	—	74 200	265	25 000
Thickwall Polyethylene*	15.0	342081	100	18 x 98	—	—	74 200	265	25 000
Thickwall Polypropylene*	15.0	342082	100	18 x 98	—	—	74 200	265	25 000
Adapters for Non-Beckman Coulter Tubes/Bottles									
Round Bottom	3.0	—	—	12 x 75	342329 (qty. of 1)	1	—	—	—

* To order caps for 15-mL tubes 342080, 342081, and 342082, use part number 343656 for a package of 50.

Rotor Replacement Parts

964911	Large O-ring for Rotor Lid
11757	Small O-ring for Rotor Lid
363024	Lid Assembly, Biosafety (Tie-down Assembly included)
363052	Rotor Lid Assembly, Single Lock (Tie-down Assembly included)

Adapters

303313 342327

